

Product Datasheet

Human IFN-gamma R1 / IFNGR1 Protein (orb257557)

Description

Human IFN-gamma R1, His Tag is expressed from



Reactivity

Gly ...

Human

Endotoxins

1.0 EU per µg

Conjugation

Unconjugated

Target

IFN-gamma R1

Preservatives

PBS, pH7.4

Form/Appearance

Powder

Storage

-20°C to -70°C for 12 months in lyophilized state; -70°C for 3 months under sterile conditions after reconstitution. For long term storage, the product should be stored at lyophilized state at -20°C or lower.

Tag

C-6×His

Note

For research use only

Application notes

This protein carries a polyhistidine tag at the C-terminus. The protein has a calculated MW of 26.6 kDa. The protein migrates as 35-45 kDa under reducing (R) condition (SDS-PAGE) due to glycosylation.

Purity

95%

Source

Human IFN-gamma R1, His Tag (orb257557) is expressed from human 293 cells (HEK293). It contains AA Glu 18 - Gly 245 (Accession # AAH05333).

MW

26.6 kDa

Uniprot ID


[P15260](#)

NCBI

[AAH05333](#)


Expiration Date

6 months from date of receipt.



SDS-PAGE analysis of Human IFN-gamma R1

...



Immobilized Human IFN-gamma, premium grade...