



Product Datasheet

Human IFN-alpha/beta R2 Protein (orb257560)

Description	Interferon alpha/beta receptor 2 (IFNAR2) is also known as IFN-alpha binding protein, IFN-alpha/beta...
Reactivity	Human
Endotoxins	1.0 EU per µg
Conjugation	Unconjugated
Target	IFN-alpha/beta R2
Preservatives	PBS, pH7.4
Form/Appearance	Powder
Storage	-20°C to -70°C for 12 months in lyophilized state; -70°C for 3 months under sterile conditions after reconstitution. For long term storage, the product should be stored at lyophilized state at -20°C or lower.
Tag	C-8×His
Note	For research use only
Application notes	This protein carries a polyhistidine tag at the C-terminus. The protein has a calculated MW of 25.8 kDa. The protein migrates as 34-45 kDa under reducing (R) condition (SDS-PAGE) due to glycosylation.
Purity	95%
Source	Human IFNAR2, His Tag (orb257560) is expressed from human 293 cells (HEK293). It contains AA Ile 27 - Lys 243 (Accession # P48551-2).
MW	25.8 kDa
NCBI	NP_997467.1
Expiration Date	6 months from date of receipt.



SDS-PAGE analysis of Human IFN-alpha/beta...



Human IFN-alpha 2b Protein, Fc Tag immob...