



Product Datasheet

Human IFN-alpha / beta R2 Protein (orb257561)

Description	Interferon alpha/beta receptor 2 (IFNAR2) is also known as IFN-alpha binding protein, IFN-alpha/beta...
Reactivity	Human
Endotoxins	1.0 EU per µg
Conjugation	Unconjugated
Target	IFN-alpha/beta R2
Preservatives	Tris with Glycine, Arginine and NaCl, pH7.5
Form/Appearance	Powder
Storage	-20°C to -70°C for 12 months in lyophilized state; -70°C for 3 months under sterile conditions after reconstitution. For long term storage, the product should be stored at lyophilized state at -20°C or lower.
Tag	C-Fc
Note	For research use only
Application notes	This protein carries a human IgG1 Fc tag at the C-terminus. The protein has a calculated MW of 51.4 kDa. The protein migrates as 66-90 kDa under reducing (R) condition (SDS-PAGE) due to glycosylation.
Purity	95%
Source	Human IFNAR2, Fc Tag (orb257561) is expressed from human 293 cells (HEK293). It contains AA Ile 27 - Lys 243 (Accession # P48551-2).
MW	51.4 kDa
NCBI	NP_997467.1
Expiration Date	6 months from date of receipt.



SDS-PAGE analysis of Human IFN-alpha/beta...



Loaded Human IFN-alpha 1, His Tag (Cat. ...