



Product Datasheet

Human IFN-gamma / IFNG Protein (orb257562)

Description	Human IFN-gamma, premium grade is expressed from human 293 cells (HEK293). It contains AA Gln 24 - G...
Reactivity	Human
Endotoxins	0.01 EU per µg
Conjugation	Unconjugated
Target	IFN-gamma
Preservatives	PBS, pH7.4
Form/Appearance	Powder
Storage	-20°C to -70°C for 12 months in lyophilized state; -70°C for 3 months under sterile conditions after reconstitution. For long term storage, the product should be stored at lyophilized state at -20°C or lower.
Tag	Native
Note	For research use only
Application notes	This protein carries no "tag". The protein has a calculated MW of 16.8 kDa. The protein migrates as 15 kDa, 19 kDa and 24 kDa±3 kDa under reducing (R) condition (SDS-PAGE) due to glycosylation.
Purity	95%
Source	Human IFN-gamma, premium grade (orb257562) is expressed from human 293 cells (HEK293). It contains AA Gln 24 - Gln 166 (Accession # AAH70256).
MW	16.8 kDa
Uniprot ID	P01579
NCBI	AAH70256
Expiration Date	6 months from date of receipt.



SDS-PAGE analysis of Human IFN gamma pro...



Immobilized Human IFN-gamma, premium gra...