



Product Datasheet

Mouse Fc Protein (orb257565)


Description	Crystallizable fragments composed of the carboxy-terminal halves of both IMMUNOGLOBULIN HEAVY CHAINS...
Reactivity	Mouse
Endotoxins	1.0 EU per µg
Conjugation	Unconjugated
Target	IgG1 Fc
Preservatives	51 mM Tris, 100 mM Glycine, 25 mM Arginine, 150 mM NaCl, pH7.5
Form/Appearance	Powder
Storage	-20°C
Tag	Native
Note	For research use only
Application notes	This protein carries no "tag". The protein has a calculated MW of 24.3 kDa. The protein migrates as 28-30 kDa under reducing (R) condition (SDS-PAGE) due to glycosylation.
Protein Sequence	AAK53870.1
Purity	95%
Source	Mouse IgG Fc, Tag Free (orb257565) is expressed from human 293 cells (HEK293). It contains AA Val 111 - Lys 324 (Accession # AAK53870.1).
MW	24.3 kDa
NCBI	AAK53870.1
Expiration Date	6 months from date of receipt.



SDS-PAGE analysis of Mouse IgG Fc protein...



Mouse IgG Fc protein



Mouse IgG Fc protein