

Product Datasheet

Human IL-1RL1 / ST2 Protein (orb257583)

Description

Human IL-1RL1, Fc Tag (IL1-H5250) is expressed from



Reactivity

Human

Endotoxins

1.0 EU per µg

Conjugation

Unconjugated

Target

IL-1RL1

Preservatives

Tris with Glycine, Arginine and NaCl, pH7.5

Form/Appearance

Powder

Storage

-20°C to -70°C for 12 months in lyophilized state; -70°C for 3 months under sterile conditions after reconstitution. For long term storage, the product should be stored at lyophilized state at -20°C or lower.

Tag

C-Fc

Note

For research use only

Application notes

This protein carries a human IgG1 Fc tag at the C-terminus. The protein has a calculated MW of 61.0 kDa. The protein migrates as 80-100 kDa under reducing (R) condition (SDS-PAGE) due to glycosylation.

Purity

● -20°C to -70°C for 12 months in lyophilized state; ● -70°C for 3 months under sterile conditions after reconstitution. For long term storage, the product should be stored at lyophilized state at -20°C or lower.

Source

Human IL-1RL1, Fc Tag (orb257583) is expressed from human 293 cells (HEK293). It contains AA Lys 19 - Pro 323 (Accession # NP_003847).

MW

61.0 kDa

Uniprot ID


[Q01638](#)

NCBI


[NP_003847](#)

Expiration Date

6 months from date of receipt.



SDS-PAGE analysis of Human IL-1RL1 prote...



Loaded Human IL-1RL1, Fc Tag (Cat. No. o...