



Product Datasheet

Mouse IL-13 R alpha 1 Protein (orb257586)

Description	Interleukin 13 receptor, alpha 1 is also known as IL13RA1, NR4 and CD213A1 (cluster of differentiati...
Reactivity	Mouse
Endotoxins	1.0 EU per µg
Conjugation	Unconjugated
Target	IL-13 R alpha 1
Preservatives	Tris with Glycine, Arginine and NaCl, pH7.5
Form/Appearance	Powder
Storage	-20°C to -70°C for 12 months in lyophilized state; -70°C for 3 months under sterile conditions after reconstitution. For long term storage, the product should be stored at lyophilized state at -20°C or lower.
Tag	C-Fc
Note	For research use only
Application notes	This protein carries a human IgG1 Fc tag at the C-terminus. The protein has a calculated MW of 62.7 kDa. The protein migrates as 67-100 kDa under reducing (R) condition (SDS-PAGE) due to glycosylation.
Purity	95%
Source	Mouse IL-13RA1, Fc Tag (orb257586) is expressed from human 293 cells (HEK293). It contains AA Ala 26 - Thr 340 (Accession # NP_598751).
MW	62.7 kDa
NCBI	NP_598751
Expiration Date	6 months from date of receipt.



SDS-PAGE analysis of Mouse IL-13 R alpha...



Loaded Mouse IL-13RA1, Fc Tag (Cat. No. ...