



Product Datasheet

Human IL-27 Ra / TCCR Protein (orb257609)

Description	Interleukin-27 receptor subunit alpha (IL27RA) is also known as Cytokine receptor WSX-1, Cytokine r...
Reactivity	Human
Endotoxins	1.0 EU per µg
Conjugation	Unconjugated
Target	IL-27 Ra
Preservatives	50 mM Tris, 100 mM Glycine, 25 mM Arginine, 150 mM NaCl, pH7.5
Form/Appearance	Powder
Storage	-20°C to -70°C for 12 months in lyophilized state; -70°C for 3 months under sterile conditions after reconstitution. For long term storage, the product should be stored at lyophilized state at -20°C or lower.
Tag	C-Fc
Note	For research use only
Application notes	This protein carries a human IgG1 Fc tag at the C-terminus. The protein has a calculated MW of 78.7 kDa. The protein migrates as 100-116 kDa under reducing (R) condition (SDS-PAGE) due to glycosylation.
Purity	95%
Source	Human IL-27 Ra, Fc Tag (orb257609) is expressed from human 293 cells (HEK293). It contains AA Gln 33 - Arg 514 (Accession # AAH28003).
MW	78.7 kDa
NCBI	AAH28003
Expiration Date	6 months from date of receipt.



SDS-PAGE analysis of Human IL-27 Ra prot...



Human IL-27 Ra Protein, Fc Tag (Cat. No....