

Product Datasheet

Smox antibody (FITC) (orb25761)

Description	Smox antibody (FITC)
Species/Host	Rabbit
Reactivity	Mouse
Conjugation	FITC
Immunogen	Recombinant Mouse Spermine oxidase protein (1-555AA)
Target	Smox
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q99K82
Dilution Range	ELISA
Expiration Date	12 months from date of receipt.