

Product Datasheet


Human IL-5 R alpha / CD125 Protein (orb257612)

Description


Human IL-5 R alpha, His Tag is expressed from

kDa	M	R

	Glu ...
Reactivity	Human
Endotoxins	1.0 EU per µg
Conjugation	Unconjugated
Target	IL-5 R alpha
Preservatives	PBS, pH7.4
Form/Appearance	Powder
Storage	-20°C to -70°C for 12 months in lyophilized state; -70°C for 3 months under sterile conditions after reconstitution. For long term storage, the product should be stored at lyophilized state at -20°C or lower.
Tag	C-10×His
Note	For research use only
Application notes	This protein carries a polyhistidine tag at the C-terminus. The protein has a calculated MW of 37.7 kDa. The protein migrates as 45-55 kDa under reducing (R) condition (SDS-PAGE) due to glycosylation.
Purity	85%
Source	Human IL-5 R alpha, His Tag (orb257612) is expressed from human 293 cells (HEK293). It contains AA Asp 21 - Glu 335 (Accession # NP_000555).
MW	37.7 kDa
Uniprot ID	Q01344
NCBI	NP_000555
Expiration Date	6 months from date of receipt.



SDS-PAGE analysis of Human IL-5 R alpha ...



Loaded Biotinylated Human IL-5, His, Avi...