



Product Datasheet

Human IL-18BP Protein (orb257623)

Description	Interleukin-18-binding protein (IL18BP) is also known as Tadekinig-alfa, IL18BPα, which functions as...
Reactivity	Human
Endotoxins	1.0 EU per µg
Conjugation	Unconjugated
Target	IL-18BP
Preservatives	Tris with Glycine, Arginine and NaCl, pH7.5
Form/Appearance	Powder
Storage	-20°C to -70°C for 12 months in lyophilized state; -70°C for 3 months under sterile conditions after reconstitution. For long term storage, the product should be stored at lyophilized state at -20°C or lower.
Tag	C-Fc
Note	For research use only
Application notes	This protein carries a human IgG1 Fc tag at the C-terminus. The protein has a calculated MW of 44.3 kDa. The protein migrates as 56-66 kDa under reducing (R) condition (SDS-PAGE) due to glycosylation.
Purity	95%
Source	Human IL-18BP, Fc Tag (orb257623) is expressed from human 293 cells (HEK293). It contains AA Thr 31 - Gly 194 (Accession # NP_001138529).
MW	44.3 kDa
NCBI	NP_001138529
Expiration Date	6 months from date of receipt.



SDS-PAGE analysis of Human IL-18BP prote...



Human IL-18BP, Fc Tag (Cat. No. orb257623)