


Product Datasheet

Human IL-4 R alpha / CD124 Protein (orb257626)

Description

Human IL-4 R alpha, His Tag is expressed from



Reactivity

His ...

Human

Endotoxins

1.0 EU per µg

Conjugation

Unconjugated

Target

IL-4 R alpha

Preservatives

PBS, pH7.4

Form/Appearance

Powder

Storage

-20°C to -70°C for 12 months in lyophilized state; -70°C for 3 months under sterile conditions after reconstitution. For long term storage, the product should be stored at lyophilized state at -20°C or lower.

Tag

C-6×His

Note

For research use only

Application notes

This protein carries a polyhistidine tag at the C-terminus. The protein has a calculated MW of 24.6 kDa. The protein migrates as 36-50 kDa when calibrated against Star Ribbon Pre-stained Protein Marker under reducing (R) condition (SDS-PAGE) due to glycosylation.

Purity

95%

Source

Human IL-4 R alpha, His Tag (orb257626) is expressed from human 293 cells (HEK293). It contains AA Met 26 - His 232 (Accession # NP_000409.1).

MW

24.6 kDa

Uniprot ID


[P24394](#)

NCBI


[NP_000409.1](#)

Expiration Date

6 months from date of receipt.



SDS-PAGE analysis of Human IL-4 R alpha ...



Immobilized Human IL-4 R alpha, His Tag ...