



Product Datasheet

Human IL-21 R Protein (orb257627)

Description	Human IL-21 R, His Tag is expressed from human 293 cells (HEK293). It contains AA Cys 20 - Pro 236 (...)
Reactivity	Human
Endotoxins	1.0 EU per µg
Conjugation	Unconjugated
Target	IL-21 R
Preservatives	PBS, pH7.4
Form/Appearance	Powder
Storage	-20°C to -70°C for 12 months in lyophilized state; -70°C for 3 months under sterile conditions after reconstitution. For long term storage, the product should be stored at lyophilized state at -20°C or lower.
Tag	C-10×His
Note	For research use only
Application notes	This protein carries a polyhistidine tag at the C-terminus. The protein has a calculated MW of 26.9 kDa. The protein migrates as 37-47 kDa under reducing (R) condition (SDS-PAGE) due to glycosylation.
Purity	90%
Source	Human IL-21 R, His Tag (orb257627) is expressed from human 293 cells (HEK293). It contains AA Cys 20 - Pro 236 (Accession # Q9HBE5-1).
MW	26.9 kDa
Uniprot ID	Q9HBE5
NCBI	NP_068570.1
Expiration Date	6 months from date of receipt.



SDS-PAGE analysis of Human IL-21 R prote...



Immobilized Biotinylated Human IL-21, Fc...