


Product Datasheet

Human LAIR2 / CD306 Protein (orb257658)

Description	Leukocyte-associated immunoglobulin-like receptor 2 (LAIR2) is also known as CD antigen CD306. LAIR2...
Reactivity	Human
Endotoxins	1.0 EU per µg
Conjugation	Unconjugated
Target	LAIR-2
Preservatives	PBS, pH7.4
Form/Appearance	Powder
Storage	-20°C
Tag	C-6×His
Note	For research use only
Application notes	This protein carries a polyhistidine tag at the C-terminus. The protein has a calculated MW of 14.9 kDa. The protein migrates as 18-20 kDa under reducing (R) condition (SDS-PAGE).
Protein Sequence	AAH69366
Purity	95%
Source	Human LAIR-2, His Tag (orb257658) is expressed from human 293 cells (HEK293). It contains AA Gln 22 - Pro 152 (Accession # AAH69366).
MW	14.9 kDa
NCBI	AAH69366
Expiration Date	6 months from date of receipt.



SDS-PAGE
analysis of
Human
LAIR-2
protei...