

Product Datasheet

Human Leptin R / CD295 Protein (orb257667)

Description Human Leptin R, Fc Tag is expressed from human 293 cells (HEK293). It contains AA Phe 22 - Asp 839 (...)

Reactivity Human

Endotoxins 1.0 EU per µg

Conjugation Unconjugated

Target Leptin R

Preservatives 50 mM Tris, 100 mM Glycine, pH7.5

Form/Appearance Powder

Storage -20°C

Tag C-Fc

Note For research use only

Application notes This protein carries a human IgG1 Fc tag at the C-terminus. The protein has a calculated MW of 119.6 kDa. The protein migrates as 150-165 kDa under reducing (R) condition (SDS-PAGE) due to glycosylation.

Protein Sequence NP_002294.2

Purity 90%


Source Human Leptin R, Fc Tag (orb257667) is expressed from human 293 cells (HEK293). It contains AA Phe 22 - Asp 839 (Accession # NP_002294.2).

MW 119.6 kDa

Uniprot ID [P48357](#)

NCBI [NP_002294.2](#)

Expiration Date 6 months from date of receipt.



SDS-PAGE
analysis of
Human
Leptin R
prot...