


Product Datasheet

Human CD58 / LFA-3 Protein (orb257668)

Description	Human CD58, His Tag is expressed from human 293 cells (HEK293). It contains AA Phe 29 - Arg 215 (Acc...
Reactivity	Human
Endotoxins	1.0 EU per µg
Conjugation	Unconjugated
Target	CD58
Preservatives	PBS, pH7.4
Form/Appearance	Powder
Storage	-20°C
Tag	C-6×His
Note	For research use only
Application notes	This protein carries a polyhistidine tag at the C-terminus. The protein has a calculated MW of 22.3 kDa. The protein migrates as 30-55 kDa under reducing (R) condition (SDS-PAGE) due to glycosylation.
Protein Sequence	AAH05930
Purity	95%
Source	Human CD58, His Tag (orb257668) is expressed from human 293 cells (HEK293). It contains AA Phe 29 - Arg 215 (Accession # AAH05930).
MW	22.3 kDa
Uniprot ID	P19256
NCBI	AAH05930
Expiration Date	6 months from date of receipt.



SDS-PAGE
analysis of
Human
CD58
protein