


Product Datasheet

Human LIF R / CD118 Protein (orb257674)


Description

Human LIF R, Fc Tag (LIR-H4252) is expressed from

Reactivity	Human
Endotoxins	1.0 EU per µg
Conjugation	Unconjugated
Target	LIF R
Preservatives	Tris with Glycine, Arginine and NaCl, pH7.5
Form/Appearance	Powder
Storage	-20°C to -70°C for 12 months in lyophilized state; -70°C for 3 months under sterile conditions after reconstitution. For long term storage, the product should be stored at lyophilized state at -20°C or lower.
Tag	C-Fc
Note	For research use only
Application notes	This protein carries a human IgG1 Fc tag at the C-terminus. The protein has a calculated MW of 115.5 kDa. The protein migrates as 140-190 kDa under reducing (R) condition (SDS-PAGE) due to glycosylation.
Purity	● -20°C to -70°C for 12 months in lyophilized state; ● -70°C for 3 months under sterile conditions after reconstitution. For long term storage, the product should be stored at lyophilized state at -20°C or lower.
Source	Human LIF R, Fc Tag (orb257674) is expressed from human 293 cells (HEK293). It contains AA Gln 45 - Ser 833 (Accession # P42702-1).
MW	115.5 kDa
Uniprot ID	P42702
NCBI	NP_001121143
Expiration Date	6 months from date of receipt.



SDS-PAGE analysis of Human LIF R protein



Immobilized Human LIF Protein, premium g...