


Product Datasheet

Human Lumican Protein (orb257678)

Description	Lumican (LUM) is also known as Keratan sulfate proteoglycan lumican, KSPG lumican, LDC and SLRR2D, wh...
Reactivity	Human
Endotoxins	1.0 EU per µg
Conjugation	Unconjugated
Target	Lumican
Preservatives	PBS, pH7.4
Form/Appearance	Powder
Storage	-20°C
Tag	C-6×His
Note	For research use only
Application notes	This protein carries a polyhistidine tag at the C-terminus. The protein has a calculated MW of 37.5 kDa. The protein migrates as 50-60 kDa under reducing (R) condition (SDS-PAGE) due to glycosylation.
Protein Sequence	NP_002336
Purity	95%
Source	Human Lumican, His Tag (orb257678) is expressed from human 293 cells (HEK293). It contains AA Gln 19 - Asn 338 (Accession # NP_002336).
MW	37.5 kDa
NCBI	NP_002336
Expiration Date	6 months from date of receipt.



SDS-PAGE
analysis of
Human
Lumican
prote...