



Product Datasheet

SNRPD1 antibody (orb25768)


Description	SNRPD1 antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IF, IHC, WB
Immunogen	Recombinant Human Small nuclear ribonucleoprotein Sm D1 protein (1-119AA)
Target	SNRPD1
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P62314
Dilution Range	WB:1:200-5000, IHC-P:1:20-200, IF:1:50-200
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293 whole ce...



Immunofluorescence analysis of PC-3 cell...



Immunohistochemical staining of human ki...