


Product Datasheet

Human LYPD3 / C4.4A Protein (orb257680)

Description

Human LYPD3, His Tag (LY3-H52H5) is expressed


Reactivity

Human

Endotoxins

1.0 EU per µg

Conjugation

Unconjugated

Target

LYPD3

Preservatives

PBS, pH7.4

Form/Appearance

Powder

Storage

-20°C

Tag

C-10×His

Note

For research use only

Application notes

This protein carries a polyhistidine tag at the C-terminus. The protein has a calculated MW of 28.7 kDa. The protein migrates as 45-66 kDa under reducing (R) condition (SDS-PAGE) due to glycosylation.

Protein Sequence

AAH39167

Purity

● -20°C to -70°C for 12 months in lyophilized state; ● -70°C for 3 months under sterile conditions after reconstitution. For long term storage, the product should be stored at lyophilized state at -20°C or lower.

Source

Human LYPD3, His Tag (orb257680) is expressed from human 293 cells (HEK293). It contains AA Leu 31 - His 286 (Accession # AAH39167).

MW

28.7 kDa

Uniprot ID
[O95274](#)
NCBI
[AAH39167](#)
Expiration Date

6 months from date of receipt.

SDS-PAGE
analysis of
Human
LYPD3
protein