


Product Datasheet

Human MICA Protein (orb257689)

Description	Human MICA, His Tag is expressed from human 293 cells (HEK293). It contains AA Glu 24 - Gln 308 (Acc...
Reactivity	Human
Endotoxins	1.0 EU per µg
Conjugation	Unconjugated
Target	MICA
Preservatives	PBS, pH7.4
Form/Appearance	Powder
Storage	-20°C
Tag	C-7×His
Note	For research use only
Application notes	This protein carries a polyhistidine tag at the C-terminus. The protein has a calculated MW of 33.8 kDa. The protein migrates as 45-66 kDa under reducing (R) condition (SDS-PAGE) due to glycosylation.
Protein Sequence	AAH16929.1
Purity	95%
Source	Human MICA, His Tag (orb257689) is expressed from human 293 cells (HEK293). It contains AA Glu 24 - Gln 308 (Accession # AAH16929.1).
MW	33.8 kDa
Uniprot ID	Q29983
NCBI	AAH16929.1
Expiration Date	6 months from date of receipt.



SDS-PAGE
analysis of
Human MICA
protein