


Product Datasheet

Human MMP-9 Protein (orb257692)

Description	Human MMP-9, His Tag is expressed from human 293 cells (HEK293). It contains AA Ala 20 - Pro 469 (Ac...
Reactivity	Human
Endotoxins	1.0 EU per µg
Conjugation	Unconjugated
Target	MMP-9
Preservatives	50 mM Tris, 150 mM NaCl, pH7.5
Form/Appearance	Powder
Storage	-20°C
Tag	C-6×His
Note	For research use only
Application notes	This protein carries a polyhistidine tag at the C-terminus. The protein has a calculated MW of 50.8 kDa. The protein migrates as 58-66 kDa under reducing (R) condition (SDS-PAGE) due to glycosylation.
Protein Sequence	AAH06093
Purity	90%
Source	Human MMP-9, His Tag (orb257692) is expressed from human 293 cells (HEK293). It contains AA Ala 20 - Pro 469 (Accession # AAH06093).
MW	50.8 kDa
Uniprot ID	P14780
NCBI	AAH06093
Expiration Date	6 months from date of receipt.



SDS-PAGE
analysis of
Human MMP-
9 protein