

Product Datasheet

Human Neprilysin / MME / CD10 Protein (orb257693)

Description Cluster of differentiation 10 (CD10) is also known as membrane metallo-endopeptidase, neutral endope...

Reactivity Human

Endotoxins 1.0 EU per µg

Conjugation Unconjugated

Target Neprilysin

Preservatives 20 mM Tris, 100 mM NaCl, pH8.0

Form/Appearance Powder

Storage -20°C

Tag Native

Note For research use only

Application notes This protein carries no "tag". The protein has a calculated MW of 80.0 kDa. The protein migrates as 80-100 kDa under reducing (R) condition (SDS-PAGE) due to glycosylation.

Protein Sequence NP_000893.2


Purity 95%

Source Human Neprilysin, Tag Free (orb257693) is expressed from human 293 cells (HEK293). It contains AA Tyr 52 - Trp 750 (Accession # NP_000893.2).

MW 80.0 kDa

NCBI [NP_000893.2](#)

Expiration Date 6 months from date of receipt.



SDS-PAGE analysis of Human Neprilysin pr...