

Product Datasheet

Human Neprilysin / MME / CD10 Protein (orb257694)

Description Human Neprilysin, Fc Tag is expressed from human 293 cells (HEK293). It contains AA Tyr 52 - Trp 750...

Reactivity Human

Endotoxins 1.0 EU per µg

Conjugation Unconjugated

Target Neprilysin

Preservatives Tris with Glycine, Arginine and NaCl, pH7.5

Form/Appearance Powder

Storage -20°C

Tag N-Fc

Note For research use only

Application notes This protein carries a human IgG1 Fc tag at the N-terminus. The protein has a calculated MW of 106.8 kDa. The protein migrates as 106-120 kDa under reducing (R) condition (SDS-PAGE) due to glycosylation.

Protein Sequence NP_000893.2

Purity 95%


Source Human Neprilysin, Fc Tag (orb257694) is expressed from human 293 cells (HEK293). It contains AA Tyr 52 - Trp 750 (Accession # NP_000893.2).

MW 106.8 kDa

Uniprot ID **P08473**

NCBI **NP_000893.2**

Expiration Date 6 months from date of receipt.



SDS-PAGE analysis of Human Neprilysin pr...