



Product Datasheet

Human NKp46 / NCR1 / CD335 Protein (orb257704)

Description	Natural cytotoxicity triggering receptor 1 (NCR1) is also known as Natural killer cell p46-related p...
Reactivity	Human
Endotoxins	0.01 EU per µg
Conjugation	Unconjugated
Target	NKp46
Preservatives	50 mM Tris, 100 mM Glycine, 25 mM Arginine, 10 mM NaCl, pH7.5
Form/Appearance	Powder
Storage	-20°C to -70°C for 12 months in lyophilized state; -70°C for 3 months under sterile conditions after reconstitution. For long term storage, the product should be stored at lyophilized state at -20°C or lower.
Tag	C-Fc
Note	For research use only
Application notes	This protein carries a human IgG1 Fc tag at the C-terminus. The protein has a calculated MW of 53.0 kDa. The protein migrates as 60-66 kDa under reducing (R) condition (SDS-PAGE) due to glycosylation.
Purity	95%
Source	Human NKp46, Fc Tag (orb257704) is expressed from human 293 cells (HEK293). It contains AA Gln 22 - Asn 254 (Accession # AAH64806).
MW	53.0 kDa
NCBI	AAH64806
Expiration Date	6 months from date of receipt.



SDS-PAGE analysis of Human NKp46 protein



Immobilized Human NKp46, Fc Tag, premium...