


Product Datasheet

Human IL-12A / NKSF1 / p35 Protein (orb257707)

Description	Human IL-12A, Fc Tag is expressed from human 293 cells (HEK293). It contains AA Arg 23 - Ser 219 (Ac...
Reactivity	Human
Endotoxins	1.0 EU per µg
Conjugation	Unconjugated
Target	IL-12A
Preservatives	Tris with Glycine, Arginine and NaCl, pH7.5
Form/Appearance	Powder
Storage	-20°C to -70°C for 12 months in lyophilized state; -70°C for 3 months under sterile conditions after reconstitution. For long term storage, the product should be stored at lyophilized state at -20°C or lower.
Tag	C-Fc
Note	For research use only
Application notes	This protein carries a human IgG1 Fc tag at the C-terminus. The protein has a calculated MW of 49.2 kDa. The protein migrates as 57-66 kDa under reducing (R) condition (SDS-PAGE) due to glycosylation.
Purity	90%
Source	Human IL-12A, Fc Tag (orb257707) is expressed from human 293 cells (HEK293). It contains AA Arg 23 - Ser 219 (Accession # P29459).
MW	49.2 kDa
Uniprot ID	P29459
Expiration Date	6 months from date of receipt.



SDS-PAGE analysis of Human IL-12A protei...