



Product Datasheet

Human IL-12B / NKSF2 / p40 Protein (orb257709)

Description	Human IL-12B, His Tag is expressed from human 293 cells (HEK293). It contains AA Ile 23 - Ser 328 (A...
Reactivity	Human
Endotoxins	1.0 EU per µg
Conjugation	Unconjugated
Target	IL-12B
Preservatives	PBS, pH7.4
Form/Appearance	Powder
Storage	-20°C to -70°C for 12 months in lyophilized state; -70°C for 3 months under sterile conditions after reconstitution. For long term storage, the product should be stored at lyophilized state at -20°C or lower.
Tag	C-10×His
Note	For research use only
Application notes	This protein carries a polyhistidine tag at the C-terminus. The protein has a calculated MW of 36.6 kDa. The protein migrates as 44 kDa and 46 kDa under reducing (R) condition (SDS-PAGE) due to glycosylation.
Purity	95%
Source	Human IL-12B, His Tag (orb257709) is expressed from human 293 cells (HEK293). It contains AA Ile 23 - Ser 328 (Accession # P29460-1).
MW	36.6 kDa
Uniprot ID	P29460
NCBI	AAH67499.1
Expiration Date	6 months from date of receipt.



SDS-PAGE analysis of Human IL-12B protei...



Loaded Human IL-12 R beta 1, Fc Tag (Cat...