


Product Datasheet


Human Neuropilin-1 / NRP1 / CD304 Protein (orb257717)

Description

Human Neuropilin-1, His Tag is expressed from



	Lys ...
Reactivity	Human
Endotoxins	1.0 EU per µg
Conjugation	Unconjugated
Target	Neuropilin-1
Preservatives	PBS, pH7.4
Form/Appearance	Powder
Storage	-20°C
Tag	C-7×His
Note	For research use only
Application notes	This protein carries a polyhistidine tag at the C-terminus. The protein has a calculated MW of 70.8 kDa. The protein migrates as 75-116 kDa under reducing (R) condition (SDS-PAGE) due to glycosylation.
Protein Sequence	AAH07533
Purity	92%
Source	Human Neuropilin-1, His Tag (orb257717) is expressed from human 293 cells (HEK293). It contains AA Phe 22 - Lys 644 (Accession # AAH07533).
MW	70.8 kDa
Uniprot ID	O14786
NCBI	AAH07533
Expiration Date	6 months from date of receipt.



SDS-PAGE
analysis of
Human
Neuropilin-
1 ...