

Product Datasheet

Human TrkB / NTRK2 Protein (orb257719)

Description

Human TrkB, His Tag (NT2-H5228) is expressed from



Reactivity Human

Endotoxins 1.0 EU per µg

Conjugation Unconjugated

Target TrkB

Preservatives PBS, pH7.4

Form/Appearance Powder

Storage -20°C

Tag C-6×His

Note For research use only

Application notes This protein carries a polyhistidine tag at the C-terminus. The protein has a calculated MW of 45.2 kDa. The protein migrates as 60-90 kDa under reducing (R) condition (SDS-PAGE) due to glycosylation.

Protein Sequence AAH31835

Purity ● -20°C to -70°C for 12 months in lyophilized state; ● -70°C for 3 months under sterile conditions after reconstitution. For long term storage, the product should be stored at lyophilized state at -20°C or lower.

Source Human TrkB, His Tag (orb257719) is expressed from human 293 cells (HEK293). It contains AA Cys 32 - His 430 (Accession # AAH31835).

MW 45.2 kDa

Uniprot ID [Q16620](#)

NCBI [AAH31835](#)

Expiration Date 6 months from date of receipt.

SDS-PAGE
analysis of
Human TrkB
protein