


Product Datasheet

Human TrkB / NTRK2 Protein (orb257720)

Description

Human TrkB, Fc Tag (NT2-H5254) is expressed from


Reactivity

Human

Endotoxins

1.0 EU per µg

Conjugation

Unconjugated

Target

TrkB

Preservatives

50 mM Tris, 100 mM Glycine, pH7.5

Form/Appearance

Powder

Storage

-20°C

Tag

C-Fc

Note

For research use only

Application notes

This protein carries a human IgG1 Fc tag at the C-terminus. The protein has a calculated MW of 70.4 kDa. The protein migrates as 86-116 kDa under reducing (R) condition (SDS-PAGE) due to glycosylation.

Protein Sequence

AAH31835

Purity

● -20°C to -70°C for 12 months in lyophilized state; ● -70°C for 3 months under sterile conditions after reconstitution. For long term storage, the product should be stored at lyophilized state at -20°C or lower.

Source

Human TrkB, Fc Tag (orb257720) is expressed from human 293 cells (HEK293). It contains AA Cys 32 - His 430 (Accession # AAH31835).

MW

70.4 kDa

Uniprot ID
[Q16620](#)
NCBI
[AAH31835](#)
Expiration Date

6 months from date of receipt.

SDS-PAGE analysis of Human TrkB protein