

Product Datasheet

Mouse OX40 / TNFRSF4 / CD134 Protein (orb257731)

Description Mouse OX40, Fc Tag is expressed from human 293 cells (HEK293). It contains AA Val 20 - Pro 211 (Acce...

Reactivity Mouse

Endotoxins 1.0 EU per µg

Conjugation Unconjugated

Target OX40

Preservatives Tris with Glycine, Arginine and NaCl, pH7.5

Form/Appearance Powder

Storage -20°C

Tag C-Fc

Note For research use only

Application notes This protein carries a human IgG1 Fc tag at the C-terminus. The protein has a calculated MW of 47.9 kDa. The protein migrates as 55-70 kDa under reducing (R) condition (SDS-PAGE) due to glycosylation.

Protein Sequence NP_035789

Purity 90%


Source Mouse OX40, Fc Tag (orb257731) is expressed from human 293 cells (HEK293). It contains AA Val 20 - Pro 211 (Accession # NP_035789).

MW 47.9 kDa

Uniprot ID [P47741](#)

NCBI [NP_035789](#)

Expiration Date 6 months from date of receipt.



SDS-PAGE
analysis of
Mouse
OX40
protein