


Product Datasheet

Human CD200 / OX-2 Protein (orb257732)

Description	CD200 is also known as OX-2 membrane glycoprotein (OX-2), is a type-1 membrane glycoprotein, which c...
Reactivity	Human
Endotoxins	1.0 EU per µg
Conjugation	Unconjugated
Target	CD200
Preservatives	PBS, pH7.4
Form/Appearance	Powder
Storage	-20°C
Tag	C-6×His
Note	For research use only
Application notes	This protein carries a polyhistidine tag at the C-terminus. The protein has a calculated MW of 23.3 kDa. The protein migrates as 35-45 kDa under reducing (R) condition (SDS-PAGE) due to glycosylation.
Protein Sequence	AAH22522
Purity	98%
Source	Human CD200, His Tag (orb257732) is expressed from human 293 cells (HEK293). It contains AA Gln 31 - Gly 232 (Accession # AAH22522).
MW	23.3 kDa
NCBI	AAH22522
Expiration Date	6 months from date of receipt.



SDS-PAGE
analysis of
Human
CD200
protein