


Product Datasheet

Human CD200 / OX-2 Protein (orb257733)

Description	Human CD200, Fc Tag (OX2-H5251) is expressed from human 293 cells (HEK293). It contains AA Gln 31 - ...
Reactivity	Human
Endotoxins	1.0 EU per µg
Conjugation	Unconjugated
Target	CD200
Preservatives	50 mM Tris, 100 mM Glycine, pH7.5
Form/Appearance	Powder
Storage	-20°C
Tag	C-Fc
Note	For research use only
Application notes	This protein carries a human IgG1 Fc tag at the C-terminus. The protein has a calculated MW of 48.6 kDa. The protein migrates as 55-68 kDa under reducing (R) condition (SDS-PAGE) due to glycosylation.
Protein Sequence	AAH22522
Purity	● -20°C to -70°C for 12 months in lyophilized state; ● -70°C for 3 months under sterile conditions after reconstitution. For long term storage, the product should be stored at lyophilized state at -20°C or lower.
Source	Human CD200, Fc Tag (orb257733) is expressed from human 293 cells (HEK293). It contains AA Gln 31 - Gly 232 (Accession # AAH22522).
MW	48.6 kDa
Uniprot ID	P41217
NCBI	AAH22522
Expiration Date	6 months from date of receipt.



SDS-PAGE
analysis of
Human
CD200
protein