



Product Datasheet

Mouse PD-L1 / B7-H1 Protein (orb257764)

Description	Mouse PD-L1, Fc Tag (PD1-M5251) is expressed from human 293 cells (HEK293). It contains AA Phe 19 - ...
Reactivity	Mouse
Endotoxins	1.0 EU per µg
Conjugation	Unconjugated
Target	PD-L1
Preservatives	Tris with Glycine, Arginine and NaCl, pH7.5
Form/Appearance	Powder
Storage	-20°C
Tag	C-Fc
Note	For research use only
Application notes	This protein carries a human IgG1 Fc tag at the C-terminus. The protein has a calculated MW of 51.4 kDa. The protein migrates as 60-80 kDa under reducing (R) condition (SDS-PAGE) due to glycosylation.
Protein Sequence	NP_068693
Purity	● -20°C to -70°C for 12 months in lyophilized state; ● -70°C for 3 months under sterile conditions after reconstitution. For long term storage, the product should be stored at lyophilized state at -20°C or lower.
Source	Mouse PD-L1, Fc Tag (orb257764) is expressed from human 293 cells (HEK293). It contains AA Phe 19 - Thr 238 (Accession # NP_068693).
MW	51.4 kDa
Uniprot ID	Q9EP73
NCBI	NP_068693
Expiration Date	6 months from date of receipt.



SDS-PAGE analysis of Mouse PD-L1 protein



Line graph illustrates about the sensitivity...