


Product Datasheet

Human PROCR / EPCR / CD201 Protein (orb257784)

Description	Endothelial protein C receptor (PROCR) is also known as Activated protein C receptor (APC receptor),...
Reactivity	Human
Endotoxins	1.0 EU per µg
Conjugation	Unconjugated
Target	PROCR
Preservatives	50 mM Tris, 100 mM Glycine, pH7.5
Form/Appearance	Powder
Storage	-20°C
Tag	C-Fc
Note	For research use only
Application notes	This protein carries a human IgG1 Fc tag at the C-terminus. The protein has a calculated MW of 48.6 kDa. The protein migrates as 60-66 kDa under reducing (R) condition (SDS-PAGE) due to glycosylation.
Protein Sequence	AAH14451
Purity	95%
Source	Human PROCR, Fc Tag (orb257784) is expressed from human 293 cells (HEK293). It contains AA Ser 18 - Thr 209 (Accession # AAH14451).
MW	48.6 kDa
NCBI	AAH14451
Expiration Date	6 months from date of receipt.



SDS-PAGE analysis of Human PROCR protein