

Product Datasheet

Human Nectin-1 / PVRL1 / CD111 Protein (orb257789)

Description

Human Nectin-1, His Tag is expressed from human


Reactivity

Human

Endotoxins

1.0 EU per µg

Conjugation

Unconjugated

Target

Nectin-1

Preservatives

PBS, pH7.4

Form/Appearance

Powder

Storage

-20°C

Tag

C-8×His

Note

For research use only

Application notes

This protein carries a polyhistidine tag at the C-terminus. The protein has a calculated MW of 35.0 kDa. The protein migrates as 45-60 kDa when calibrated against Star Ribbon Pre-stained Protein Marker under reducing (R) condition (SDS-PAGE) due to glycosylation.

Protein Sequence

NP_002846

Purity

95%

Source

Human Nectin-1, His Tag (orb257789) is expressed from human 293 cells (HEK293). It contains AA Gln 31 - Thr 334 (Accession # NP_002846).

MW

35.0 kDa

Uniprot ID

[Q15223](#)

NCBI

[NP_002846](#)

Expiration Date

6 months from date of receipt.

SDS-PAGE
analysis of
Human
Nectin-1
prot...