


Product Datasheet

Mouse RENIN Protein (orb257794)

Description	Mouse RENIN, His Tag is expressed from human 293 cells (HEK293). It contains AA Leu 22 - Arg 402 (Ac...
Reactivity	Mouse
Endotoxins	1.0 EU per µg
Conjugation	Unconjugated
Target	RENIN
Preservatives	PBS, pH7.4
Form/Appearance	Powder
Storage	-20°C
Tag	C-6×His
Note	For research use only
Application notes	This protein carries a polyhistidine tag at the C-terminus. The protein has a calculated MW of 37.0 kDa. The protein migrates as 47-57 kDa under reducing (R) condition (SDS-PAGE) due to glycosylation.
Protein Sequence	NP_112469.1
Purity	95%
Source	Mouse RENIN, His Tag (orb257794) is expressed from human 293 cells (HEK293). It contains AA Leu 22 - Arg 402 (Accession # NP_112469.1).
MW	37.0 kDa
Uniprot ID	P00797
NCBI	NP_112469.1
Expiration Date	6 months from date of receipt.



SDS-PAGE
analysis of
Mouse RENIN
protein