


Product Datasheet

Mouse ROR1 Protein (orb257797)

Description	Tyrosine-protein kinase transmembrane receptor ROR1 is also known as Neurotrophic tyrosine kinase, r...
Reactivity	Mouse
Endotoxins	0.1 EU per µg
Conjugation	Unconjugated
Target	ROR1
Preservatives	PBS, pH7.4
Form/Appearance	Powder
Storage	-20°C
Tag	C-8×His
Note	For research use only
Application notes	This protein carries a polyhistidine tag at the C-terminus. The protein has a calculated MW of 42.9 kDa. The protein migrates as 55-65 kDa under reducing (R) condition (SDS-PAGE) due to glycosylation.
Protein Sequence	Gln 30 - Glu 403
Purity	95%
Source	Mouse ROR1, His Tag (orb257797) is expressed from human 293 cells (HEK293). It contains AA Gln 30 - Glu 403 (Accession # Q9Z139).
MW	42.9 kDa
Expiration Date	6 months from date of receipt.



SDS-PAGE analysis of Mouse ROR1 protein