


Product Datasheet

Mouse ROR1 Protein (orb257798)

Description	Tyrosine-protein kinase transmembrane receptor ROR1 is also known as Neurotrophic tyrosine kinase, r...
Reactivity	Mouse
Endotoxins	1.0 EU per µg
Conjugation	Unconjugated
Target	ROR1
Preservatives	50 mM Tris, 100 mM Glycine, 25 mM Arginine, 150 mM NaCl, pH7.5
Form/Appearance	Powder
Storage	-20°C
Tag	C-Fc
Note	For research use only
Application notes	This protein carries a human IgG1 Fc tag at the C-terminus. The protein has a calculated MW of 68.4 kDa. The protein migrates as 76-96 kDa under reducing (R) condition (SDS-PAGE) due to glycosylation.
Protein Sequence	Gln 30 - Glu 403
Purity	95%
Source	Mouse ROR1, Fc Tag (orb257798) is expressed from human 293 cells (HEK293). It contains AA Gln 30 - Glu 403 (Accession # Q9Z139).
MW	68.4 kDa
Expiration Date	6 months from date of receipt.



SDS-PAGE analysis of Mouse ROR1 protein