



Product Datasheet

Human ROR2 / NTRKR2 Protein (orb257799)

Description	Tyrosine-protein kinase transmembrane receptor ROR2 is also known as Neurotrophic tyrosine kinase, r...
Reactivity	Human
Endotoxins	1.0 EU per µg
Conjugation	Unconjugated
Target	ROR2
Preservatives	Tris with Glycine, Arginine and NaCl, pH7.5
Form/Appearance	Powder
Storage	-20°C
Tag	C-Fc
Note	For research use only
Application notes	This protein carries a human IgG1 Fc tag at the C-terminus. The protein has a calculated MW of 68 kDa. The protein migrates as 70-90 kDa under reducing (R) condition (SDS-PAGE) due to glycosylation.
Protein Sequence	AAI30523.1
Purity	95%
Source	Human ROR2, Fc Tag (orb257799) is expressed from human 293 cells (HEK293). It contains AA Glu 34 - Gly 403 (Accession # A1L4F5-1).
MW	68 kDa
NCBI	AAI30523.1
Expiration Date	6 months from date of receipt.



SDS-PAGE analysis of Human ROR2 protein



Line graph illustrates about the Ag-Ab r...