


Product Datasheet

Human R-Spondin 1 / RSPO1 Protein (orb257807)

Description	R-spondin-1 is also known as Roof plate-specific Spondin 1 (RSPO1) and cysteine-rich and single throm...
Reactivity	Human
Endotoxins	1.0 EU per µg
Conjugation	Unconjugated
Target	R-Spondin 1
Preservatives	PBS, pH7.4
Form/Appearance	Powder
Storage	-20°C
Tag	C-10×His
Note	For research use only
Application notes	This protein carries a polyhistidine tag at the C-terminus. The protein has a calculated MW of 28.7 kDa. The protein migrates as 38-42 kDa under reducing (R) condition (SDS-PAGE) due to glycosylation.
Protein Sequence	NP_001033722.1
Purity	95%
Source	Human R-Spondin 1, His Tag (orb257807) is expressed from human 293 cells (HEK293). It contains AA Ser 21 - Ala 263 (Accession # Q2MKA7-1).
MW	28.7 kDa
NCBI	NP_001033722.1
Expiration Date	6 months from date of receipt.



SDS-PAGE
analysis of
Human R-
Spondin 1
p...