


Product Datasheet

Human S100A6 Protein (orb257817)

Description	S100 calcium binding protein A6 (S100-A6) is also known as 2A9, 5B10, CABP, CACY and PRA, is a member...
Reactivity	Human
Endotoxins	1.0 EU per µg
Conjugation	Unconjugated
Target	S100A6
Preservatives	PBS, pH7.4
Form/Appearance	Powder
Storage	-20°C
Tag	Native
Note	For research use only
Application notes	This protein carries no "tag". The protein has a calculated MW of 10.2 kDa. The protein migrates as 10.2 kDa under reducing (R) condition (SDS-PAGE).
Protein Sequence	AAH01431
Purity	95%
Source	Human S100A6, Tag Free (orb257817) is expressed from E.coli cells. It contains AA Met 1 - Gly 90 (Accession # AAH01431).
MW	10.2 kDa
NCBI	AAH01431
Expiration Date	6 months from date of receipt.



SDS-PAGE
analysis of
Human
S100A6
protei...