


Product Datasheet

Human S100A8 Protein (orb257818)

Description	S100 calcium binding protein A8 (S100-A8) is also known as Cystic fibrosis antigen, Migration inhibit...
Reactivity	Human
Endotoxins	1.0 EU per µg
Conjugation	Unconjugated
Target	S100A8
Preservatives	PBS, pH7.4
Form/Appearance	Powder
Storage	-20°C
Tag	Native
Note	For research use only
Application notes	This protein carries no "tag". The protein has a calculated MW of 10.8 kDa. The protein migrates as 12 kDa under reducing (R) condition (SDS-PAGE).
Protein Sequence	AAH05928
Purity	95%
Source	Human S100A8, Tag Free (orb257818) is expressed from E.coli cells. It contains AA Met 1 - Glu 93 (Accession # AAH05928).
MW	10.8 kDa
NCBI	AAH05928
Expiration Date	6 months from date of receipt.



SDS-PAGE
analysis of
Human
S100A8
protei...