


Product Datasheet

Human S100P / S100E Protein (orb257820)

Description	Human S100P, Tag Free is expressed from E. coli cells. It contains AA Met 1 - Lys 95 (Accession # AA...
Reactivity	Human
Endotoxins	1.0 EU per µg
Conjugation	Unconjugated
Target	S100P
Preservatives	50 mM Tris, 150 mM NaCl, pH8.0
Form/Appearance	Powder
Storage	-20°C
Tag	Native
Note	For research use only
Application notes	This protein carries no "tag". The protein has a calculated MW of 10.4 kDa. The protein migrates as 10.4 kDa under reducing (R) condition (SDS-PAGE).
Protein Sequence	AAH06819
Purity	95%
Source	Human S100P, Tag Free (orb257820) is expressed from E.coli cells. It contains AA Met 1 - Lys 95 (Accession # AAH06819).
MW	10.4 kDa
Uniprot ID	P25815
NCBI	AAH06819
Expiration Date	6 months from date of receipt.



SDS-PAGE
analysis of
Human
S100P
protein