

Product Datasheet

Human Serpin F1 / PEDF Protein (orb257828)

Description Human Serpin F1, His Tag is expressed from human 293 cells (HEK293). It contains AA Gln 20 - Pro 418...

Reactivity Human

Endotoxins 1.0 EU per µg

Conjugation Unconjugated

Target Serpin F1

Preservatives PBS, pH7.4

Form/Appearance Powder

Storage -20°C

Tag C-6×His

Note For research use only

Application notes This protein carries a polyhistidine tag at the C-terminus. The protein has a calculated MW of 45.2 kDa. The protein migrates as 45-60 kDa under reducing (R) condition (SDS-PAGE) due to glycosylation.

Protein Sequence NP_002606

Purity 95%


Source Human Serpin F1, His Tag (orb257828) is expressed from human 293 cells (HEK293). It contains AA Gln 20 - Pro 418 (Accession # NP_002606).

MW 45.2 kDa

Uniprot ID **P36955**

NCBI **NP_002606**

Expiration Date 6 months from date of receipt.



SDS-PAGE
analysis of
Human
PEDF
protein