


Product Datasheet

Human Serpin A3 Protein (orb257833)

Description	Serpin A3 is also known as Alpha-1-antichymotrypsin (ACT or AACT), Cell growth-inhibiting gene 24/25...
Reactivity	Human
Endotoxins	1.0 EU per µg
Conjugation	Unconjugated
Target	Serpin A3
Preservatives	PBS, pH7.4
Form/Appearance	Powder
Storage	-20°C
Tag	C-10×His
Note	For research use only
Application notes	This protein carries a polyhistidine tag at the C-terminus. The protein has a calculated MW of 46.6 kDa. The protein migrates as 50-63 kDa under reducing (R) condition (SDS-PAGE) due to glycosylation.
Protein Sequence	AAH03559
Purity	95%
Source	Human Serpin A3, His Tag (orb257833) is expressed from human 293 cells (HEK293). It contains AA His 24 - Ala 423 (Accession # AAH03559).
MW	46.6 kDa
Uniprot ID	P01011
NCBI	AAH03559
Expiration Date	6 months from date of receipt.



SDS-PAGE
analysis of
Human
Serpin A3
pro...