

## Product Datasheet

### Human Semaphorin 4A / SEMA4A Protein (orb257835)

## Description

Human Semaphorin 4A, His Tag is expressed from



## Reactivity

His...

Human

## Endotoxins

1.0 EU per µg

## Conjugation

Unconjugated

## Target

Semaphorin 4A

## Preservatives

PBS, pH7.4

## Form/Appearance

Powder

## Storage

-20°C

## Tag

C-6×His

## Note

For research use only

## Application notes

This protein carries a polyhistidine tag at the C-terminus. The protein has a calculated MW of 72.7 kDa. The protein migrates as 70-80 kDa under reducing (R) condition (SDS-PAGE) due to glycosylation.

## Protein Sequence

AAH20974

## Purity

95%

## Source

Human Semaphorin 4A, His Tag (orb257835) is expressed from human 293 cells (HEK293). It contains AA Gly 33 - His 683 (Accession # AAH20974).

## MW

72.7 kDa

## Uniprot ID


[Q9H3S1](#)

## NCBI

[AAH20974](#)

## Expiration Date

6 months from date of receipt.



SDS-PAGE  
analysis of  
Human  
Semaphorin  
4A...