

Product Datasheet

Human Siglec-5 / CD170 Protein (orb257840)

Description Human Siglec-5, Fc Tag is expressed from human 293 cells (HEK293). It contains AA Glu 17 - Thr 434 (...)

Reactivity Human

Endotoxins 1.0 EU per µg

Conjugation Unconjugated

Target Siglec-5

Preservatives Tris with Glycine, Arginine and NaCl, pH7.5

Form/Appearance Powder

Storage -20°C

Tag C-Fc

Note For research use only

Application notes This protein carries a human IgG1 Fc tag at the C-terminus. The protein has a calculated MW of 73.0 kDa. The protein migrates as 90-110 kDa under reducing (R) condition (SDS-PAGE) due to glycosylation.

Protein Sequence AAH29896

Purity 95%


Source Human Siglec-5, Fc Tag (orb257840) is expressed from human 293 cells (HEK293). It contains AA Glu 17 - Thr 434 (Accession # AAH29896).

MW 73.0 kDa

Uniprot ID [O15389](#)

NCBI [AAH29896](#)

Expiration Date 6 months from date of receipt.



SDS-PAGE
analysis of
Human
Siglec-5
prot...