

Product Datasheet

Human Siglec-6 / CD327 Protein (orb257841)

Description Human Siglec-6, Fc Tag is expressed from human 293 cells (HEK293). It contains AA Gln 16 - Val 320 (...)

Reactivity Human

Endotoxins 1.0 EU per µg

Conjugation Unconjugated

Target Siglec-6

Preservatives 50 mM Tris, 100 mM Glycine, 150 mM NaCl, pH7.5

Form/Appearance Powder

Storage -20°C

Tag C-Fc

Note For research use only

Application notes This protein carries a human IgG1 Fc tag at the C-terminus. The protein has a calculated MW of 60.3 kDa. The protein migrates as 80-100 kDa under reducing (R) condition, and 130-150 kDa and 250-300 kDa under non-reducing (NR) condition (SDS-PAGE) due to glycosylation.

Protein Sequence AAH35359

Purity 95%


Source Human Siglec-6, Fc Tag (orb257841) is expressed from human 293 cells (HEK293). It contains AA Gln 16 - Val 320 (Accession # AAH35359).

MW 60.3 kDa

Uniprot ID [O43699](#)

NCBI [AAH35359](#)

Expiration Date 6 months from date of receipt.



SDS-PAGE
analysis of
Human
Siglec-6
prot...