

Product Datasheet

Human SIRP alpha / CD172a Protein (orb257844)

Description Human SIRP alpha, Mouse IgG1 Fc Tag is expressed from human 293 cells (HEK293). It contains AA Glu 3...

Reactivity Human

Endotoxins 1.0 EU per µg

Conjugation Unconjugated

Target SIRP alpha

Preservatives 50 mM Tris, 100 mM Glycine, pH7.5

Form/Appearance Powder

Storage -20°C

Tag C-mFc(mIgG1)

Note For research use only

Application notes This protein carries a mouse IgG1 Fc tag at the C-terminus. The protein has a calculated MW of 63.0 kDa. The protein migrates as 70-105 kDa under reducing (R) condition (SDS-PAGE) due to glycosylation.

Protein Sequence NP_001035111

Purity 95%


Source Human SIRP alpha, Mouse IgG1 Fc Tag (orb257844) is expressed from human 293 cells (HEK293). It contains AA Glu 31 - Arg 370 (Accession # NP_001035111).

MW 63.0 kDa


Uniprot ID **P78324**

NCBI **NP_001035111**

Expiration Date 6 months from date of receipt.



SDS-PAGE analysis of Human SIRP alpha pr...



Line graph illustrates about the sensiti...